

COMPANY INTRODUCTION

SHAOXING UNDIS MEDICAL TECHNOLOGY CO., LTD is a leading manufacturer of good-quality disposable anesthesia products.

Everything for Health. This is not only our company's slogan, but also our attitude for controlling the quality of products. All our devices are manufactured in compliance with ISO9001, ISO13485 and National Medical Industry Standard. With CE & FDA marking, we could sell our goods freely throughout European and the world with no restriction. We always regard quality as our life by insisting on independent innovation and scientific research. We have our R&D center, advanced productive and testing equipment and is also equipped with class 100,000 clean workshop.

Products is our capital. We are committed to manufacturing a wide range of innovative and high-quality anesthesia disposables, such as Anesthesia breathing circuit(corrugated circuit, expandable circuit, smooth bore circuit), breathing filter(bacterial filter, HMEF Filter), Catheter Mount, Anesthesia Mask, Laryngeal Mask Airway, Breathing Bag, etc. And we are continuing to develop new products.

Customers' trust is our wealth. Our goal is to continually supply high-quality anesthesia disposable for hospitals and clinics, meeting the market demand both at home and abroad. Through our unremitting efforts, our products have been exported to many countries, such as France, USA, Mexico, Peru, Turkey, Iran, Vietnam, India, Libya, Thailand, etc. Till now, we have earned good reputation.



“Everything For Health”




Disposable Corrugated Breathing Circuits

Features:

1. Corrugated tubes have a good flexibility and resistance to bend.
2. Consist of corrugated hoses, swivel elbow with luer port and cap, Y-piece connector.
3. According to ISO standard, match all kinds of anesthetic and breathing machines.
4. ISO, CE and FDA certifications.




1111602	Adult basic corrugated parallel 1.6m	20	
1121602	Pediatric basic corrugated parallel 1.6m	20	



1111612	Adult basic corrugated parallel 1.6m with 0.9m extra limb	20	
1121613	Pediatric basic corrugated parallel 1.6m with 0.9m extra limb	20	



1111605	Adult basic corrugated parallel 1.6m with 0.9m extra limb and 2L latex free breathing bag	20	
1121605	Pediatric basic corrugated parallel 1.6m with 0.9m extra limb and 1L latex free breathing bag	20	





1111604	Adult basic corrugated parallel 1.6m with two water traps	20	
1121604	Pediatric basic corrugated parallel 1.6m with two water traps	20	




Disposable Smooth Bore Breathing Circuits

Features:


1. Smooth bore breathing systems incorporate flexible end connectors at the machine end and light weight rigid connectors within the system offering improved ventilatory support.
2. The smooth bore tube can improve flow characteristics and lower tubing compliance.
3. ISO, CE and FDA certifications.




1311602	Adult smooth bore 1.6m	20	
1321602	Pediatric smooth bore 1.6m	20	



1311605	Adult smooth bore 1.6m with 0.9m extra limb	20	
1321605	Pediatric smooth bore 1.6m with 0.9m extra limb	20	



1311612	Adult smooth bore 1.6m with 0.9m extra limb and 2L latex free breathing bag	20	
1321613	Pediatric smooth bore 1.6m with 0.9m extra limb and 1L latex free breathing bag	20	




1311604	Adult smooth bore 1.6m with two water traps	20	
1321604	Pediatric smooth bore 1.6m with two water traps	20	

Disposable Expandable Breathing Circuits

Features:

1. Adjustable length can significantly save space.
2. Expandable tube can be packed more in per box, will lower freight volume costs.
3. Customized tube colors and lengths are accepted.
4. ISO, CE and FDA certifications.




1211602	Adult expandable circuit 1.6m	20	
1221602	Pediatric expandable circuit 1.6m	20	



1211612	Adult expandable circuit 1.6m with 0.9m extra limb	20	
1221613	Pediatric expandable circuit 1.6m with 0.9m extra limb	20	



1211605	Adult expandable circuit 1.6m with 0.9m extra limb and 2L latex free breathing bag	20	
1221605	Pediatric expandable circuit 1.6m with 0.9m extra limb and 1L latex free breathing bag	20	




1211604	Adult expandable circuit 1.6m with two water traps	20	
1221604	Pediatric expandable circuit 1.6m with two water traps	20	

Disposable Coaxial Circuits

Features:


1. Coaxial breathing circuit is a convenient "tube within a tube" breathing system designed with low flow anaesthesia in mind, offering warmer and more humid inspiratory gases, and the option of an integrated respiratory gases monitoring line.
2. Available in a variety of sizes:
3. According to ISO Standard, match all kinds of anaesthetic and breathing machines.
4. Customized tube colors and lengths are accepted.
5. ISO, CE and FDA certifications



1411603

Coaxial circuit 1.6m with 0.4m expandable tube

20



1411616

Coaxial circuit 1.6m with 0.9m extra limb and 2L latex free breathing bag

20

Disposable Duo-limb Circuits

Features:


1. Less weight than two-limb circuits, reduces torque on the patient's airway.
2. With one single limb, offers more versatility when used as transport circuit and in the OR.
3. According to ISO Standard, match all kinds of anaesthetic and breathing machines.
4. Customized tube colors and lengths are accepted.
5. ISO, CE and FDA certifications



1511602

Duo-limb circuit 1.6m with 0.4m expandable tube

20



1521605

Duo-limb circuit 1.6m with 0.9m extra limb and 2L latex free breathing bag

20



Heat & Moisture Exchange Filter (HMEF)

High Filtration Efficiency help to protect against Bacterial & Virus contamination and by humidifying the inhaled air, patient' s airways can be protected from cold and dry gases.

Model	3000400	Certificate	CE / ISO
Dimensions	Height:77mm width:73mm	Dead space	20ml
Weight	28g	Tidal Volume	200-1500ml
Moisture Output	25mg@500VT	Connector	22M/15F-22F/15M
Humidity loss	15mg	Materials	Shell: PP Membrane: HV Brand, Electrostatic
Filtration Efficiency	BFE 99.9999% VFE 99.9999%	Pressure Drop @ 30LPM	Dry: 120pa Wet: 150pa
Package	Paper-plastic bag	Units	100pcs
Carton size	38-32-34cm	Weight	3.5kg/2.8kg



Model	3000500	Certificate	CE / ISO
Dimensions	Height:79mm width:64mm	Dead space	20ml
Weight	22g	Tidal Volume	200-1500ml
Moisture Output	25mg@500VT	Connector	22F-22M/15F
Humidity loss	15mg	Materials	Shell: PP Membrane: HV Brand, Electrostatic
Filtration Efficiency	BFE 99.9999% VFE 99.9999%	Pressure Drop @ 30LPM	Dry: 120pa Wet: 150pa
Package	Paper-plastic bag	Units	100pcs
Carton size	38-32-34cm	Weight	3.5kg/2.8kg



Model	3000600	Certificate	CE / ISO
Dimensions	Height:66mm width:40mm	Dead space	12ml
Weight	13.5g	Tidal Volume	150-300ml
Moisture Output	24mg@500VT	Connector	22F-22M/15F
Humidity loss	12mg	Materials	Shell: PP Membrane: HV Brand, Electrostatic
Filtration Efficiency	BFE 99.9999% VFE 99.9999%	Pressure Drop @ 30LPM	Dry: 120pa Wet: 140pa
Package	Paper-plastic bag	Units	100pcs
Carton size	38-32-34cm	Weight	3.5kg/2.8kg



Model	3000700	Certificate	CE / ISO
Dimensions	Height:51mm width:38mm	Dead space	8ml
Weight	8g	Tidal Volume	70-150ml
Moisture Output	22mg@500VT	Connector	15F-15M/8M
Humidity loss	10mg	Materials	Shell: PP Membrane: HV Brand, Electrostatic
Filtration Efficiency	BFE 99.9999% VFE 99.9999%	Pressure Drop @ 30LPM	Dry: 120pa Wet: 140pa
Package	Paper-plastic bag	Units	100pcs
Carton size	38-32-34cm	Weight	2.0kg/1.3kg





Bacterial & Virus Filter

Connected with breathing circuit, help patients to prevent from Cross-contamination by providing 99.999% filtration efficiency of bacteria and virus.

Model	3000100	Certificate	CE / ISO
Dimensions	Height:74mm width:67mm	Dead space	25ml
Weight	22g	Tidal Volume	200-1500ml
Materials	Shell: PP Membrane: HV Brand,Electrostatic	Connector	22F-22M/15F
Filtration Efficiency	BFE 99.9999% VFE 99.9999%	Pressure Drop @ 30LPM	Dry: 120pa Wet: 140pa
Package	Paper-plastic bag	Units	100pcs
Carton size	38-32-34cm	Weight	3.5kg/2.8kg



Model	3000200	Certificate	CE / ISO
Dimensions	Height:74mm width:67mm	Dead space	25ml
Weight	22g	Tidal Volume	200-1500ml
Materials	Shell: PP Membrane: HV Brand,Electrostatic	Connector	22F-22M/15F
Filtration Efficiency	BFE 99.9999% VFE 99.9999%	Pressure Drop @ 30LPM	Dry: 120pa Wet: 140pa
Package	Paper-plastic bag	Units	100pcs
Carton size	38-32-34cm	Weight	3.5kg/2.8kg



Model	3000300	Certificate	CE / ISO
Dimensions	Height:87mm width:74mm	Dead space	35ml
Weight	53g	Tidal Volume	200-1500ml
Materials	Shell: PP Membrane: HV Brand,Electrostatic	Connector	22F-22M/15F
Filtration Efficiency	BFE 99.9999% VFE 99.9999%	Pressure Drop @ 30LPM	Dry: 80pa Wet: 100pa
Package	Paper-plastic bag	Units	100pcs
Carton size	50-45-37cm	Weight	6.1kg/5.3kg





Heat & Moisture Exchanger

As a part of Breathing system in Operation and ICU, Maximize the Humidity and Heat for the inhaled gas and protect the patient from the harm of cold and dry air.

Model	3000800	Certificate	CE / ISO
Dimensions	Height:28mm width:36mm	Pressure Drop @ 30LPM	50 Pa
Weight	4.5g	Tidal Volume	200-1500ml
Materials	Shell: PP Membrane:Sponge	Connector	15F
Moisture Loss	15mg	Moisture Output	24mg@500VT
Package	Paper-plastic bag	Units	300pcs
Carton size	63-36-48cm	Weight	4.2/3.2kg



Model	3000900	Certificate	CE / ISO
Dimensions	Height:46mm width:34mm	Pressure Drop @ 30LPM	50 Pa
Weight	9.5g	Tidal Volume	200-1500ml
Materials	Shell: pp Medium: HME Paper	Connector	15F
Moisture Loss	18mg	Moisture Output	25mg@500VT
Package	Paper-plastic bag	Units	300pcs
Carton size	63-36-48cm	Weight	5.7/4.7kg



Model	3001000	Certificate	CE / ISO
Dimensions	Height:28mm width:38mm	Pressure Drop @ 30LPM	50 Pa
Weight	5g	Tidal Volume	200-1500ml
Materials	Shell: pp Medium: HME Paper	Connector	15F
Moisture Loss	15mg	Moisture Output	25mg@500VT
Package	Paper-plastic bag	Units	300pcs
Carton size	63-36-48cm	Weight	4.2/3.2kg





Catheter Mounts

A series of catheter mounts with different tubes, elbows, ports, caps and connectors.


Expandable Range



2215110 Φ 15mm Expandable Tube + Double Swivel Elbow and Elastomeric Cap + 15M Connector + Red Cap 100



2215111 Φ 15mm Expandable Tube + Double Swivel Elbow and Elastomeric Cap + 22M/15F Connector + Red Cap 100



2215112 Φ 15mm Expandable Tube + Double Swivel Elbow and Elastomeric Cap + 22F Connector + Red Cap 100



2215113 Φ 15mm Expandable Tube + Detachable Double Swivel Elbow and Elastomeric Cap + 22M/15F Connector + 22F Connector 100



2215120 Φ 15mm Expandable Tube + Detachable Low Profile Elbow + 22M/15F Connector + 22F Connector 100



2215121 Φ 15mm Expandable Tube + Detachable Luer Lock Elbow with Cap + 22M/15F Connector + 22F Connector 100



2215122 Φ 15mm Expandable Tube + Low Profile Elbow + 22F Connector 100



2215123 Φ 15mm Expandable Tube + Luer Lock Elbow with Cap + 22F Connector 100



2215130 Φ 15mm Expandable Tube + 15M connector + 22M/15F connector + Plug 100



2215131 Φ 15mm Expandable Tube + 22M/15F Connector + 22M/15F Connector 100



2215140 Φ 15mm Expandable Tube + Neonate HMEF + 15M-15M Connector + 15M-22M/15F Connector + Plug 100



2215141 Φ 15mm Expandable Tube + Pediatric HMEF + 15M-15M Connector + 15M-22M/15F Connector + Plug 100



2215142 Φ 15mm Expandable Tube + Adult HMEF + 15M-15M Connector + 15M-22M/15F Connector + Plug 100

Smooth Bore Range



2315110 Φ 15mm Smoothbore Tube + Double Swivel Elbow and Elastomeric Cap + 15M Connector + Red Cap 100





2315111 Φ 15mm Smoothbore Tube + Double Swivel Elbow and Elastomeric Cap + 22F Connector + Red Cap 100

Smooth Bore Range

	
2315112 Φ 15mm Smoothbore Tube + Double Swivel Elbow and Elastomeric Cap + 22M/15F Connector + Red Cap 100	2315113 Φ 15mm Smoothbore Tube + Detachable Double Swivel Elbow and Elastomeric Cap + 22M/15F Connector + 22F Connector 100
	
2315120 Φ 15mm Smoothbore Tube + Detachable Luer Lock Elbow with Cap + 22M/15F Connector + 22F Connector 100	2315121 Φ 15mm Smoothbore Tube + Detachable Low Profile Elbow + 22M/15F Connector + 22F Connector 100
	
2315122 Φ 15mm Smoothbore Tube + Luer Lock Elbow with Cap + 22F Connector 100	2315123 Φ 15mm Smoothbore Tube + Low Profile Elbow + 22F Connector 100
	
2315130 Φ 15mm Smoothbore Tube + 15M connector + 22M/15F connector + Plug 100	2315131 Φ 15mm Smoothbore Tube + 22M/15F Connector + 22F connector 100

Corrugated Range

	
2115110 Φ 15mm Corrugated Tube + Double Swivel Elbow and Elastomeric Cap + 15M Connector + Red Cap 100	2115111 Φ 15mm Corrugated Tube + Double Swivel Elbow and Elastomeric Cap + 22M/15F Connector + Red Cap 100
	
2115112 Φ 15mm Corrugated Tube + Double Swivel Elbow and Elastomeric Cap + 22F Connector + Red Cap 100	2115113 Φ 15mm Corrugated Tube + Detachable Double Swivel Elbow and Elastomeric Cap + 22M/15F Connector + 22F Connector 100
	
2115120 Φ 15mm Corrugated Tube + Detachable Luer Lock Elbow with Cap + 22M/15F Connector + 22F Connector 100	2115121 Φ 15mm Corrugated Tube + Detachable Low Profile Elbow 22M/15F Connector + 22F Connector 100
	
2115122 Φ 15mm Corrugated Tube + Low Profile Elbow + 22F Connector 100	2115123 Φ 15mm Corrugated Tube + Luer Lock Elbow with Cap + 22F Connector 100
	
2115130 Φ 15mm Corrugated Tube + 15M connector + 22M/15F connector + Plug 100	2115131 Φ 15mm Corrugated Tube + 22M/15F Connector + 22F connector 100

Three-way Connector



6000101
15F/22M-22MX2



6000102
15F/22M-22MX2



6000103
15F/22M-22MX2



6000104
15F/22M-22MX2



6000105
15F/22M-22MX2



6000106
15F/22M-22MX2



6000107
22M/15F-22F/15F

Elbow Connector



6000201
22M/15F-15M



6000202
22M/15F-15M



6000203
22M/15F-15M



6000204
22M/15F-15M

Straight Connector



6000301
22M-22F



6000302
15M-22F



6000303
22F



6000304
22F

Straight Connector



6000305
22M-22M/15F



6000306
22M-22M



6000307
15M-22M/15F



6000308
15M-22F



6000309
15M-15F



6000310
15M-15M



6000311
15M-15M

Safety Cap



6000401
Luer Port Cap



6000402
Luer Port Cap



6000403
22F Connector Cap



6000404
Safety Cap

Accessories



6000501
22F




6000502
22F



6000503
22M



7000110 Disposable Anesthesia Mask
Size: 0, 1, 2, 3, 4, 5, 6




7000120 Disposable Anesthesia Mask
Size: 0, 1, 2, 3, 4, 5, 6




7000130 Disposable Nebulizer Mask
Size: S,M,L,XL



7000210 Disposable Laryngeal Mask Airway-PVC
Size: 1#,1.5#,2#,2.5#,3#,4#,5#



7000220 Disposable Laryngeal Mask Airway-Silicon
Size: 1#,1.5#,2#,2.5#,3#,4#,5#



7000310 Disposable Endotracheal Tube
Size: 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0,6.5, 7.0, 7.5, 8.0, 8.5, 9.0



7000320 Disposable Endotracheal Tube-Reinforce
Size: 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0



7000410 Disposable Oropharyngeal Airway/Guedel Airway
Size: 4,5,6,7,8,9,10,11



7000510 Disposable Nasopharyngeal Airway
Size: 6.0,6.5,7.0,7.5,8.0,8.5



7000610 Disposable Tracheotomy Tube
Size: 3.0, 3.5, 4.0, 4.5.....9.5, 10.0